


# OraLine®

## Single Step Oral Saliva Drug Test



OraLine® detects the presence of the parent THC compound rather than the metabolized drug. OraLine® is able to detect very low levels of THC over a longer period of time than is possible with competitive tests.



### Why OraLine and not the competition?

As the graph illustrates, other oral fluid drug screens on the market are simply not sensitive enough to detect marijuana for 24 hours following use.

Drug Class	Cutoff Level (ng/ml)
THC – Cannabinoids delta 9 (Marijuana, Hash, Pot)	4ng/ml
COC – Cocaine (Crack, Coke)	25ng/ml
Benzoylcegonine (COC metabolite)	10ng/ml
MET – Methamphetamine (Speed, Crank, Ecstasy)	50ng/ml
OPI – Opiates (Heroin, Morphine)	40ng/ml